


Short ACQ Introduction


- experimental setup
 - detectors, electronics, data acquisition
- event structure
 - file format, data types
- sorting into ROOT files
 - batch files, spectra, calibration

Experimental Setup




- CACTUS
- SiRi preamplifiers
- CACTUS high voltage
- SiRi high voltage

Electronics




- CACTUS amplifiers
- SiRi amplifiers
- trigger logic
- ADCs and TDCs
- VME computer

Trigger Logic






DATA Acquisition

- VME computer (bobcat):
 - eventbuilder collects data from ADCs, ...
- Linux PC (tiger):
 - `usb-engine` receives data from VME computer and writes to disk
 - `acq_master` controls engine
 - `rupdate` shows spectra
 - `mesytec-rc.py` adjusts amplifier settings



Event Structure



... ..



- each buffer is 128kB large
- event header: C0 00 00 nn
- data:
 - b = box identifier, 7 bit
 - c = channel number, 7 bit
 - v = value, 16 bit
- next event's header
- ...
- end of buffer: 80 00 00 00

Box and Channel Numbers

- TPU, E, ΔE should always be present
- Nal **should** have t and E at the same time
- clock and scaler
 - every $\sim 1000/4000$ events
 - in pairs of 2 channels

TPU	0x00	0
clock	0x01	16, 17
scaler	0x02	0..31
Nal time 0..27	0x10	0..31
Nal energy 0..27	0x24	0..31
E	0x21	0..31
ΔE 0..3	0x22	0..31
ΔE 4..7	0x23	0..31

Event Data (1)

- TPU: which trigger – presently unused information
- E back detector chips:
 - guard rings: 0, 2, ..., 14
 - detectors: 1, 3, ... 15
- ΔE : front detectors
 - 0..7 for back detector 0
 - 8..15 for back detector 1
 - ...

Event Data (2)

- Nal time:
 - time from “back over threshold” until “Nal over threshold”
 - channels 0..27
- Nal energy:
 - channels 0..27

Other Data

- clock:
 - seconds since 1.1.1970 in 2×16 bit
 - every ~ 1024 events
- scalers:
 - ΔE_{1+2} : scaler 0 = channels 0+1
(high 16 bit + low 16 bit)
 - NaI 0..15: scaler 1 = channels 2+3
 - NaI 16..27: scaler 2 = channels 4+5
 - E: scaler 3 = channels 6+7
 - every ~ 4096 events

Sorting-Program for ROOT

- **copy** `userroot.h`, `userroot.cpp`, `Makefile`, `user_routine_root.cpp`
 - **user directory** in `.../sirius/src/sort`
- **modify** `Makefile` and, if necessary, `user_routine_root.cpp`
- **run** `make` to compile `offline_sort`
- **run** `offline_sort XY.batch` to sort
- **look at spectra** in ROOT

Batch Files for Sorting

- # ... = comment
- gain file ... = load gainshift from file
- parameter ... = set parameter value
- outfile ... = set output filename
- data file ... = read event file


```
# batch file for 92Zr before TDC range change
gain file gainshifts.dat
parameter tnai_corr_enai = 12.985 30447.6 -58.8935 -4.43951e-03
parameter tnai_corr_esi = 12.641 -1.08751e+05 -516.198 -9.95697e-05
parameter ex_from_ede = 16551.9 -0.969197 -9.57725e-07 \
    16540.8 -0.968965 -9.66113e-07 \
    16528.9 -0.968702 -9.74982e-07 \
    16515.9 -0.968367 -9.8629e-07 \
    16502 -0.967961 -1.00006e-06 \
    16487.1 -0.96749 -1.01574e-06 \
    16471.1 -0.966949 -1.03386e-06 \
    16453.9 -0.966319 -1.05469e-06
#parameter ex_corr_exp = 0 1 0 1 0 1 0 1 0 1 0 1 0 1
parameter ex_corr_exp = -116.2561 1.0124 \
-122.8256 1.0118 \
-120.0651 1.0097 \
-120.0082 1.0078 \
-118.2018 1.0051 \
-114.3768 1.0019 \
-109.4958 0.9984 \
-104.1663 0.9945
outfile tmp/sorted_92zr_pp.root
data file datafiles/sirius-20090607-131232.data.gz
...
```

ROOT Spectra

- E: ΔE for all strips:
 - `m_e_de_b0f0, ..., m_e_de_b7f7, ...`
- NaI time:Nai energy
 - individual: `m_nai_e_t_01, ..., m_nai_e_t_28`
 - all summed up: `m_nai_e_t`
 - corrected: `m_nai_e_t_c`
- Nai time:E back
 - for each back dectector: `m_siri_e_t_0, ..., m_siri_e_t_7`
 - all together: `m_siri_e_t`


SiRi Calibration

- use ^{12}C data with sufficient statistics
- calculate energies for ground state and first excited state with rkinz
- find all $3 \times 8 \times 8$ peak positions in E and DE
 - “only” 384 peaks
 - ROOT macro peaks2D.C might help
- extract calibration coefficients by fitting
 - spreadsheet might help



Nal Calibration



- again, use ^{12}C data with good statistics
- find 511keV and 4.4MeV peaks for all 28 Nal detectors in time:energy matrices (m_nai_e_t_00, ...)
 - ROOT macro peaks2D.c might help
- extract calibration coefficients by fitting
 - spreadsheet might help



Nal Time Correction

- determine points on curve in `m_nai_e_t`
 - ROOT: toolbar, new Graph, click, click, click, ...


```
{ for(int i=0; i<Gname->GetN(); ++i) { cout << Gname->GetX()[i] << ' ' << Gname->GetY()[i] << endl; } }
```
- fit curve with function:
 - `TF1* f3=new TF1("f3", "200+[0]+[1]/(x+[2])+[3]*x")`
 - this goes into the parameter `tnai_corr_enai`
- same for E_{back} with `m_siri_e_t`
 - parameter `tnai_corr_esi`



Other Steps

- check that everything is okay
 - e.g., by looking at m...evol where various things are plotted as a function of time
- determine excitation energy
 - usually need re-fit to get correct energies
- make ALFNA matrix
 - use macro to write ROOT matrix in MAMA format
- run mama
- ...