



KJM 9250

**$^1\text{H}$ - $^{15}\text{N}$  f3 Experiments on the AVI and AVII-600 Spectrometers**

Version 4.0

Topspin 1.3 Windows XP AVI600

Topspin 2.1 Windows 7 AVII600



© Professor Emeritus Alistair Lawrence Wilkins,  
University of Waikato, New Zealand.  
January 2020

# **<sup>1</sup>H-<sup>15</sup>N f3 Experiments on the AVI and AVII-600 Spectrometers**

## **1.1 Introduction**

aw coded **<sup>1</sup>H-<sup>15</sup>N f3 HSQC, HSQC-DIPS12 and HMBC** parameter sets on the AVI and AVII-600 spectrometers are set up with 2048 (2K) acquired and processed points acquired across a **12 ppm <sup>1</sup>H** window centered at **4.7 ppm** and a **30 ppm <sup>15</sup>N window (SW)** centered at **118 ppm (O3)**. The <sup>1</sup>H and <sup>15</sup>N NMR signals of peptides and microcystins typically occur in these windows. Different **SW** and **O3** settings may (will) be required for other nitrogen containing compounds.

**D24 = 1/8J** (for all NH<sub>x</sub>'s, x = 1 or 2) and **D26 = 1/4J** are auto calculated from **CNST4 = <sup>1</sup>J <sup>1</sup>H-<sup>15</sup>N coupling constant = 90 Hz** in aw coded **f3 HSQC** parameter sets.

## **1.2 Processing**

**f3 HSQC** and **HSQC-DIPS12** experiments are phase sensitive experiments which should be phased *before* using the **abs1** and **abs2** commands.

**f3 HMBC** experiments are absolute value experiment. Phasing is not required.

## **2.0 Experiments and Parameter Sets**

The following aw coded **<sup>1</sup>H-<sup>15</sup>N f3 HSQC, HSQC-DIPS12 and HMBC** parameter sets are available on the **AVI-600** and AVII-600 spectrometers.

- 2.1 <sup>1</sup>H-<sup>15</sup>N f3 HSQC** spectrum
- 2.2 <sup>1</sup>H-<sup>15</sup>N f3 HSQC-DIPS12** spectrum
- 2.3 <sup>1</sup>H-<sup>15</sup>N f3 HMBCET** spectrum

## 2.1 $^1\text{H}$ - $^{15}\text{N}$ f3 HSQC spectrum

Parameter set: **awf3hsqc (+ getprosol)**

Pulse programme: **awf3hsqc**

Type **eda** (enter) and review the following default parameters:

**SW  $^1\text{H}$**  = 12 ppm, **SW  $^{15}\text{N}$**  = 30 ppm (or other suitable values).

**TD  $^1\text{H}$**  = 2K, **TD  $^{15}\text{N}$**  = 128-160 (your choice).

**O1P** =  $^1\text{H}$  spectral window midpoint = 4.7 ppm other value of your choice.

**O3P** =  $^{15}\text{N}$  spectral window midpoint = 118 ppm other value of your choice.

**NS** = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

**D1** = repetition delay = 1.0 sec or other time of your choice.

**CNST4** =  $^1J$   $^{15}\text{N}$ - $^1\text{H}$  coupling constant = 90 Hz or other value of your choice.

**D24** =  $1/8J$  (for all  $\text{NH}_x$ 's) and **D26** =  $1/4J$  are auto calculated from **CNST4**

**ZGOPTNS** = Not used.

Check gradient settings are OK for  $^{15}\text{N}$ .


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, abs1 and abs2**



600 MHz  $^1\text{H}$ - $^{15}\text{N}$  f3 HSQC spectrum of a peptide that has 9 amino acid units. 600 MHz

## 2.2 $^1\text{H}$ - $^{15}\text{N}$ f3 HSQC-DIPSII2 spectrum

Parameter set: **awf3hsqc-dipsi2 (+ getprosol)**

Pulse programme: **awf3hsqdi2f3gpsi**

Type **eda** (enter) and review the following default parameters:

**SW  $^1\text{H}$**  = 12 ppm, **SW  $^{15}\text{N}$**  = 30 ppm (or other suitable values).

**TD  $^1\text{H}$**  = 2K, **TD  $^{15}\text{N}$**  = 128-160 (your choice).

**O1P** =  $^1\text{H}$  spectral window midpoint = 4.7 ppm other value of your choice.

**O3P** =  $^{15}\text{N}$  spectral window midpoint = 118 ppm other value of your choice.

**NS** = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters including:

**D1** = repetition delay = 1.0 sec or other time of your choice.

**CNST4** =  $^1\text{J } ^{15}\text{N}$ - $^1\text{H}$  coupling constant = 90 Hz or other value of your choice.

**D24** =  $1/8\text{J}$  (for all  $\text{NH}_x$ 's) and **D26** =  $1/4\text{J}$  are auto calculated from **CNST4**


**D9** = 80 msec or other value of your choice (20-160 msec).

**ZGOPTNS** = Not used.

Check gradient settings are OK for  $^{15}\text{N}$ .

Set receiver gain using **RGA** (Important!).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**  
**WDW(F1) = WDW(F2) = QSINE**  
**SSB(F2) = SSB(F1) = 2**  
**xfb, abs1 and abs2**



600 MHz  $^1\text{H}$ - $^{15}\text{N}$  f3 HSQC-DIPSII2 spectrum of a peptide that has 9 amino acid units. The spectrum was acquired with  $d_9 = 120$  msec (rather than 80 msec).

### 2.3 $^1\text{H}$ - $^{15}\text{N}$ f3 HMBCEt spectrum

Parameter set: **awf3hmbcet+ getprosol)**

Pulse programme: **hmbcetf3gpnd**

Type **eda** (enter) and review the following default parameters

**SW  $^1\text{H}$**  = 12 ppm, **SW  $^{15}\text{N}$**  = **30** ppm (or other suitable values).

**TD  $^1\text{H}$**  = 2K, **TD  $^{15}\text{N}$**  = 96-160 (your choice).

**O1P** =  $^1\text{H}$  spectral window midpoint = 4.7 ppm other value of your choice.

**O3P** =  $^{15}\text{N}$  spectral window midpoint = 118 ppm other value of your choice.

**NS** = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

**D1** = repetition delay = **1.0 sec** or other time of your choice.

**CNST4** =  $^1\text{J}$   $^{15}\text{N}$ - $^1\text{H}$  coupling constant = **90 Hz** or other value of your choice.

**CNST13** = **6 Hz**

**ZGOPTNS** = Not used.

Check gradient settings are OK for  $^{15}\text{N}$ .


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m + abs1 and abs2**



600 MHz  $^1\text{H}$ - $^{15}\text{N}$  f3 HMBCEt spectrum of a peptide that has 9 amino acid units

### 3.0 Appendix - Some Processing Options

#### 3.1 Qfil mode processing

The vertical axis noise pattern often seen in 2D spectra when spectra are run in **9:1 H<sub>2</sub>O:D<sub>2</sub>O** or other solvents can be suppressed by setting up the experiment with its **O1 Hz**, or **O1P ppm** value (typically ~ 4.7 ppm), to that of the **HOD** line and processing the spectrum using the following **ProcPars (edp)** settings:

**BC\_MOD = qfil**

**BCFW = 0.5 ppm** or other (smaller) suppression band width value of your choice.

**COROFFS (Hz)** can be used to offset the center of the concealed region from O1

Default values are: **BC\_MOD = quad** or **no**, **BCFW = 0** or **1.00000**, **COROFFS = 0**

Baseline correction			
ABSG	5	5	Degree of polynomial for abs (0..5)
ABSF1 [ppm]	1000.00000	1000.00000	Left limit for absf
ABSF2 [ppm]	-1000.00000	-1000.00000	Right limit for absf, abs1, abs2
BCFW [ppm]	0.50000	1.00000	Filter width for bc (sfil/qfil)
COROFFS [Hz]	0	0	Correction offset for BC_MOD=spol etc.
BC_mod	qfil	no	Fid baseline modes for em, ft, xfb,...

#### 3.2 Linear prediction and STSI processing

Provided the s/n ratio of Fourier transformed <sup>15</sup>N axis data points is reasonable **linear prediction** can be applied to improve the resolution of correlations in that axis.

Fourier transform			
TDefF	0	0	Number of fid data points used by ft
STSR	0	0	First output point of strip transform
STSI	1024	0	Total number of output points of strip transform
ME_mod	no	LPfc	Linear prediction for ft, xfb, ...
NCOEF	0	32	Number of LP coefficients
LPBIN	0	256	Number of output points for LP
TDoff	0	0	Number of back-predicted points

**F1 axis (2nd column) <sup>15</sup>N settings** can be set up as:

**ME\_MOD = LPfc**, **NCOEF = 32**, **LPBIN =** twice the number of acquired increments.

Default linear prediction values when they are not used are:

**ME\_MOD = no**, **NCOEF = 0**, **LPBIN = 0**

**STSI** can be used to ONLY display spectral data to the left hand side (higher ppm side) of a selected number of processed points. If, for example, a noisy residual H<sub>2</sub>O/HOD line appeared in the vicinity of **4.7 ppm** in a **2048 point** processed spectrum acquired with **SW = 12 ppm** and **O1 = 6 ppm**, it would not be visible in the **6-12 ppm** region view of the spectrum that would be displayed when it was processed with **STSI = 1024 points**.

The default value of **STSI** is **0**.