



KJM 9250

AVII-600 ^{13}C NMR spectra.

Version 7.3

Topspin 3.2 Windows 7 AVII600



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2020

AVII-600 ¹³C NMR spectra

1.0 Introduction

All 600 MHz **aw coded** carbon parameter files have 64 K points across a 240 ppm window with D1 typically = 1.5 or 4 sec for quaternary carbon only spectra.

*Quaternary carbons may have moderate to long T1's. If saturation is suspected **D1** should be increased from its default time.*

Shaped pulse DEPT, INEPT, JMOD and APT spectra have wider flat topped spectra windows. This is most noticeable for signals within 0-20 ppm ppm of the high and low field sides of 240 ppm window ¹³C spectra, especially so in INEPT, JMOD and APT spectra.

When processed with **EF** or **EFP** a default **LB** of **1 or 2 Hz** will be applied.

2.0 ¹³C NMR Spectra

- 2.1 ¹³C spectra with power gated ¹H decoupling and NOE
- 2.2 Inverse gated ¹³C spectra (no NOE)
- 2.3 Coupled ¹³C NMR spectrum with NOE
- 2.4 ZRESTSE, ZRESTSEIG and ZRESTSEND spectra
- 2.5 DEPT45, DEPT90 and DEPT135 spectra
DEPT45SP, DEPT90SP and DEPT135SP spectra
¹H Coupled DEPTND spectra
- 2.6 DEPTQ135 spectrum with quaternary carbons
DEPTQ quaternary (singlet) carbons only spectrum
- 2.7 INEPT spectra
INEPTSP spectra
INEPTND spectrum
- 2.8 JMOD spectra
JMODSP spectra
- 2.9 APT spectra
APTSP spectra

2.1 ^{13}C NMR spectra using a 30°, 45°, 70° or 90° pulse

Parameter sets: **awcarbon30**, **awcarbon45**, **awcarbon70** or **awcarbon90** (+ **getprosol**)
Pulse programmes: **zgpg30/awzgp30**, **awzgp45**, **awzgp70** or **zgpg/awzgp90** respectively.
Spectra are ^1H decoupled with NOE.


TD = 64 K, **SI** = 64 K.
SW = 240 ppm, **O1P** = 110 ppm.
NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec or other time of your choice.
DE = 50 or 60 usec will reduce but not eliminate baseline roll in cyroprobe ^{13}C spectra.

Type ased (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies **LB**). A large negative **PHC1** phase correction (eg -650°) will be required. After phasing baseline roll can be eliminated using the **multiabsn** macro with n = 30-40.



Lower: ^{13}C spectrum of quinine in D_6 -DMSO with a 70° pulse.

Upper: ^{13}C spectrum of quinine in D_6 -DMSO with a 30° pulse.

2.2 Inverse gated ^{13}C NMR spectra

Parameter sets: **awcarbon30ig**, **awcarbon45ig**, **awcarbon70ig**, **awcarbon90ig** (+ **getprosol**)

Pulse programmes: **awzgif30**, **awzgif45**, **awzgif70**, **awzgif** respectively

Spectra are ^1H decoupled without NOE

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.


D1 = 1.5 sec or other time of your choice.

DE = 50 or 60 usec will reduce but not eliminate baseline roll in cyroprobe ^{13}C spectra.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies **LB**). A large negative **PHC1** phase correction (eg -650°) will be required. After phasing baseline roll can be eliminated using the **multiabsn** macro with n = 30-50.



Lower: Inverse gated ^{13}C spectrum of quinine in D_6 -DMSO with a 70° pulse.

Upper: Inverse gated ^{13}C spectrum of quinine in D_6 -DMSO with a 30° pulse.

2.3 ^1H coupled ^{13}C NMR spectrum with NOE

Parameter set: **awcarbon70nd (+ getprosol)**

Pulse programme: **awzg70nd**

^1H decoupled ^{13}C spectrum with NOE using a 70° pulse

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.


D1 = 1.5 sec or other time of your choice.

DE = 50 or 60 usec will reduce but not eliminate baseline roll in cyroprobe ^{13}C spectra.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies **LB**). A large negative **PHC1** phase correction (eg -650°) will be required. After phasing baseline roll can be eliminated using the **multiabsn** macro with $n = 30-50$.



Lower: ^{13}C spectrum of quinine in $\text{D}_6\text{-DMSO}$ with a 70° pulse.

Upper: ^1H coupled ^{13}C spectrum of quinine in $\text{D}_6\text{-DMSO}$ with a 70° pulse.

2.4 ZRESTSE ¹³C spectra

Parameter sets: **awzrestse**, **awzrestseig**, **awzrestsend** (+ **getprosol**)

Pulse programmes: **awzrestse**, **awzrestseig**, **awzrestsend** respectively

Spectra are free from baseline roll.

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

P0 = **P1*****CNST0/90.0**; **CNST0** = **70°** or other tip angle of your choice (30°-90° range).


D1 = **1.5 sec** (zrestse and zrestseig) or **4 sec** (zrestsend) or other time of your choice.

Type **ased** (enter) and review parameters used in the job.

Check shaped pulses and gradients are OK.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies **LB**).



ZRESTSE spectra of quinine in D₆-DMSO determined with a 70° P0 pulse.

Lower: power gated **zrestse** spectrum. **Center:** inverse gated **zrestseig** spectrum.

Upper: ¹H coupled **zrestsend** spectrum.

2.5.1 DEPT45, DEPT90 and DEPT135 spectra

Parameter sets: **awdept45**, **awdept90** or **awdept135** (+ **getprosol**)

Pulse programmes: **dept45**, **dept90** or **dept135** respectively

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.


D1 = 1.5 sec or other time of your choice.

CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **AFP** (applies LB).



DEPT spectra of quinine in D₆-DMSO. **Lower: DEPT45** spectrum.

Center: DEPT90 spectrum. **Upper: DEPT135** spectrum.

2.5.2 DEPTSP45, DEPTSP90 and DEPTSP135 spectra

Parameter sets: **awdept45sp**, **awdept90sp** or **awdept135sp** (+ **getprosol**)

Pulse programmes: **deptsp45**, **deptsp90** or **deptsp135** respectively

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.


D1 = 1.5 sec or other time of your choice.

CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).



DEPTSP spectra of quinine in D₆-DMSO. **Lower: DEPT45SP** spectrum.

Center: DEPT90SP spectrum. **Upper: DEPT135SP** spectrum.

2.5.3 ^1H Coupled DEPTND Spectra

Parameter set: **awdeptnd (+ getprosol)**

Pulse programme: **deptnd**

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec or other time of your choice.


CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST12 = 0.5, 1.0 or 1.5 for a 45, 90 or 135 degree pulse

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **AFP** (applies LB).



Lower: DEPT45 spectrum of quinine in D_6 -DMSO.

Upper: DEPT45ND spectrum.

2.6.1 DEPT135Q spectrum with quaternary carbons

Parameter set: **awdept135q (+ getprosol)**

Pulse programme: **deptqgppsp**

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec or other time of your choice.

CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.


CNST12 = **1.5** for DEPT135 spectrum with quaternary carbons.

Type **ased** (enter) and review parameters used in the job.

Check **GPZ1**, **GPZ2** and **GPZ3** gradients = **31%**.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).



Lower: ^{13}C NMR spectra of quinine in D_6 -DMSO. **Upper:** Dept135 spectrum.
CH and CH_3 carbons positive; C (quaternary) and CH_2 carbons negative.

2.6.2 DEPT Quaternary carbon only spectrum

Parameter set: **awdeptq (+ getprosol)**

Pulse programme: **deptqgppsp**

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = multiple of 4 or 8, **DS** = 8 or 16.

D1 = 1.5 sec or other time of your choice.

CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.


CNST12 = 1.5 for alternating cancellation of DEPT135 carbon signals.

Type **ased** (enter) and review parameters used in the job.

Check **GPZ1**, **GPZ2** and **GPZ3** gradients = **31%**, **31%** and **11%** respectively.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).



Lower: ¹³C NMR spectra of quinine in D₆-DMSO. **Upper:** DeptQ spectrum.

2.7.1 INEPT45, INEPT90 or INEPT135 spectra

Parameter sets: **awinept45**, **awinept90**, **awinept135** (+ **getprosol**)
Pulse programme: **ineptird** with **CNST11 = 6, 4 or 3** respectively

TD = 64 K, **SI** = 64 K.
SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec or other time of your choice.


CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST11 = **6** for INEPT45,
 = **4** for INEPT90
 = **3** for INEPT135

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).



INEPT spectra of quinine in D₆-DMSO. **Lower:** INEPT45 spectrum.
Center: INEPT90 spectrum. **Upper:** INEPT135 spectrum.

2.7.2 INEPT45SP, INEPT90SP or INEPT135SP spectra

Parameter sets: **awinept45sp**, **awinept90sp**, **awinept135sp** (+ **getprosol**)

Pulse programme: **ineptdsp** with **CNST11 = 6, 4 or 3** respectively

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec or other time of your choice.

CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST11 = **6** for INEPT45SP,


= **4** for INEPT90SP

= **3** for INEPT135SP

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB)



INEPTSP spectra of quinine in D_6 -DMSO. **Lower:** INEPT145SP spectrum.

Center: INEPT90SP spectrum **Upper:** INEPT135SP spectrum.

2.7.3 INEPTND spectrum

Parameter set: **awineptnd (+ getprosol)**

Pulse programme: **ineptnd**

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec or other time of your choice.

CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.


Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).

Signals have non-binomial intensities: CH (d) = (-1, 1); CH₂ (t) = (-1, 0, 1);

CH₃ (q) = (-1, -1, 1, 1)



INEPTND spectrum of quinone in D₆-DMSO. The multiplicity of some signals is shown.

Insert: Expansions of the 20-30 ppm regions of **ineptnd**, **zrestsend** and **inept135sp** spectra with signal annotations.

2.8.1 JMOD (J-modulated) spectra

Parameter sets: **awjmod** or **awjmodq** (+ **getprosol**)

Pulse programme: **jmod** with **cnst11 = 1** or **2**

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec (jmod) or 4 sec (jmodq) spectra or other time of your choice.


CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST11 = **1** (CH, CH₃ positive, C, CH₂ negative), or
= **2** for C (quaternary) only carbons.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*)

Process with **EF** or **EFP** (applies LB).



JMOD (lower) and **JMODQ (upper)** spectra of quinine in D₆-DMSO. Some residual positive or negatively phased signals are seen ex protonated carbons whose 1J coupling constants differ significantly from **CNST2** = 145 Hz.

2.8.2 JMODSP spectra with shaped refocusing pulses

Parameter set: awjmodsp or awjmodqsp (+ getprosol)

Pulse programme: **awjmodsp** with **cnst 11 = 1 or 2** respectively

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = multiple of 4 or 8, **DS** = 4, 8 or 16.

D1 = 1.5 sec (jmodsp) or 3-4 sec (jmodqsp) spectra or other time of your choice.


CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST11 = **1** (CH, CH₃ positive, C, CH₂ negative), or
= **2** for C (quaternary) only carbons

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).



JMODSP (lower) and **JMODSPQ (upper)** spectra of quinine in D₆-DMSO. Some residual positive or negatively phased signals are seen ex protonated carbons whose 1J coupling constants differ significantly from **CNST2** = 145 Hz.

2.9.1 APT (Attached Proton Test) spectra

Parameter sets: **awapt** or **awaptq** (+ **getprosol**)

Pulse programme: **apt** with **cnst 11 = 1** or **2** respectively

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = Any number, **DS** = 4, 8 or 16.

D1 = 1.5 sec (APT) or 3-4 sec (APTQ) spectra or other time of your choice.


CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST11 = **1** (CH, CH₃ positive, C, CH₂ negative), or
= **2** for C (quaternary) only carbons.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important*).

Process with **EF** or **EFP** (applies LB).



APT (lower) and **APTQ (upper)** spectra of quinine in D₆-DMSO. Some residual positive or negatively phased signals are seen ex protonated carbons whose 1J coupling constants differ significantly from **CNST2** = 145 Hz.

2.9.2 APTSP spectra with shaped refocusing pulses

Parameter sets: **awaptsp** or **awaptqsp** (+ **getprosol**)

Pulse programme: **awaptsp** with **cnst11** = **1** or **2** respectively

TD = 64 K, **SI** = 64 K.

SW = 240 ppm, **O1P** = 110 ppm.

NS = multiple of 4 or 8, **DS** = 4, 8 or 16.

P0 = 90° or other tip angle of your choice (30° - 90°).

D1 = 1.5 sec (ATPSP) or 4 sec (APTQSP) spectra or other time of your choice.


CNST2 = 145 Hz = $^1J(^{13}\text{C}-^1\text{H})$ or other value of your choice.

CNST11 = **1** (CH, CH₃ positive, C, CH₂ negative), or
= **2** for C (quaternary) only carbons.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EF** or **EFP** (applies LB).



APTSP (lower) and **APTSPQ (upper)** spectra of quinine in D₆-DMSO. Some residual positive or negatively phased signals are seen ex protonated carbons whose 1J coupling constants differ significantly from **CNST2** = 145 Hz.