



KJM 9250

AVII-600 SHSQC and SHMBC Experiments

Version 7.3

Topspin 3.2 Windows 7 AVII600



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2020

AVII-600 SHSQC and SHMBC Experiments

1.0 Introduction

aw coded **AVII-600 SHSQC** and **SHMBC** parameter sets are set up with 2K acquired ^1H points and 128 ^{13}C increments. ^1H and ^{13}C spectral windows and their mid points should be determined before setting up **SHSQC** or **SHMBC** experiments.

The ^{13}C axis resolution of **SHSQC** and **SHMBC** spectra acquired using 128 increments is typically 4-5 times greater than that of standard full window **HSQC** and **HMBC** spectra acquired with 160-256 or more increments.

A greater range of **SHSQC** and **SHMBC** experiments is available on the **AVII-600** than is the case for the case for the **AVI-600**

1.2 Processing

SHSQC45 and **SHSQC135** experiments are phase sensitive experiments. These spectra should be phased **before** using the **abs1** and **abs2** commands. Low level 2J correlations may be observed in **SHSQC** spectra.

SHMBCQ5 spectra are processed with **xfb**. **SHMBCCT** spectra are acquired in echo-antiecho mode and transformed with **xfb** *and* **xf2m**.

2.0 AVII-600 Experiments and Parameter Sets

2.1	SHSQC45	semi-selective hsqc45
2.2	SHSQC45.m	medium-selective hsqc45
2.3	SHSQC135	semi-selective hsqc135
2.4	SHSQC135.m	medium-selective hsqc135
2.5	SHMBCQ5	semi-selective shmbc experiment
2.6	SHMBCCT	constant time mode shmbc experiment

2.1 SHSQC45

Parameter set: **awshsqc45 (+ getprosol)**

Pulse programme: **awshsqcetgpsisp2.2-45**

d24 is automatically calculated from **cnst2**

Type **eda** (enter) and enter **SW(¹H)** and **O1P = ¹H** spectral window midpoint in ppm.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

SW(¹³C) = 20 ppm, excitation band width is ca **10-12 ppm**.

TD(F2) = 2K, **TD(F1) = 128** points.

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **3000 usec p43:sp32** band selective **Q3.1000** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024** points

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



AVII-600 SHSQC45 spectrum of quinine in D₆-DSMO. The ¹³C axis is centered at 4 ppm.

2.2 SHSQC45.m

Parameter set: **awshsqc45.m (+ getprosol)**

Pulse programme: **awshsqcetgpsisp2.2-45m**

d24 is automatically calculated from **cnst2**

Type **eda** (enter) and enter **SW(¹H)** and **O1P = ¹H** spectral window midpoint in ppm.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

SW(¹³C) = 60 ppm, excitation band width is ca **40-45 ppm**.

TD(F2) = 2K, **TD(F1) = 128** points.

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **700 usec p33:sp23** band selective **Q3.1000** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024** points

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



AVII-600 SHSQC45.m spectrum of quinine in D₆-DSMO. The ¹³C axis is centered at 40 ppm.

2.3 SHSQC135

Parameter set: **awshsqc135 (+ getprosol)**

Pulse programme: **awshsqcedetgpsisp2.3-135**

d21 and **d24** are automatically calculated from **cnst2**

Type **eda** (enter) and enter **SW(¹H)** and **O1P = ¹H** spectral window midpoint in ppm.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

SW(¹³C) = 20 ppm, excitation band width is ca **10-12 ppm**

TD(F2) = 2K, **TD(F1) = 64-128** points.

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **3000 usec p24:sp32** band selective **Q3.1000** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024** points

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



AVII-600 SHSQC135 spectrum of quinine in D₆-DSMO. The ¹³C axis is centered at 40 ppm.

2.4 SHSQC135.m

Parameter set: **awshsqc135.m (+ getprosol)**

Pulse programme: **awshsqcedetgpsisp2.3-135m**

d21 and **d24** are automatically calculated from **cnst2**

Type **eda** (enter) and enter **SW(¹H)** and **O1P = ¹H** spectral window midpoint in ppm.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

SW(¹³C) = 60 ppm, excitation band width is ca **40-45 ppm**

TD(F2) = 2K, **TD(F1) = 64-128** points.

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **768 usec p30:sp23** band selective **Q3.1000** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024** points

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



AVII-600 SHSQC135.m spectrum of quinine in D₆-DSMO. The ¹³C axis is centered at 40 ppm.

2.5 SHMBCQ5

Parameter set: **awshmbcq5 (+ getprosol)**

Pulse programme: **awshmbcq5**

Type **eda** (enter) and enter **SW(¹H)** and **O1P = ¹H** spectral window midpoint in ppm.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

SW(¹³C) = 20 ppm, excitation band width is ca **12 ppm**

TD(F2) = 2K, **TD(F1) = 64-128** points

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg 125-220 Hz).

CNST13 = ⁿJ selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **3000 usec p35:sp27** band selective **Q5.1000** pulse.


Set receiver gain using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024** points

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



AVII-600 SHMBCQ5 spectrum of quinine in D₆-DMSO. The ¹³C axis is centered at 122 ppm.

2.6 SHMBCCT (Constant time mode SHMBC experiment)

Parameter set: **awshmbcct (+ getprosol)**

Pulse programme: **awshmbcctetgpl2nd**

Type **eda** (enter) and enter **SW(¹H)** and **O1P = ¹H** spectral window midpoint in ppm.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

SW(¹³C) = 20 ppm, excitation band width is ca **10-12 ppm**.

TD(F2) = 2K, **TD(F1) = 64-128** points

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST6 = 120 Hz, **CNST7 = 170 Hz** = min/max ¹J coupling constants.

CNST13 = ⁿJ selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **3000 usec p43:sp32** band selective **Q3.1000** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024** points

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb *and* **xf2m**, **abs1** and **abs2**



AVII-600 SHMBCCT spectrum of quinine in D₆-DMSO. The ¹³C axis is centered at 122 ppm.