



KJM 9250

AVII-600 SELHSQC and SELHMBC Experiments

Version 7.3

Topspin 3.2 Windows 7 AVII600



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2020

AVII-600 SELHSQC and SELHMBC Experiments

1.0 Introduction

The ^{13}C **O1 frequency** of a target signal must be determined before running a **SELHSQC** or **SELHMBC** experiment and entered as the **O2 frequency** when setting up these experiments.

SELHSQC spectra should be processed as per standard ^1H NMR spectra.

SELHMBCQ5 spectra should be processed with **EFP *and* MC**.

2.0 SELHSQC and SELHMBC Experiments

2.1 SELHSQC and SELHSQCND spectra

2.2 SELHSQC-DIPS12 and SELHSQCND-DIPS12 spectra

2.3 SELHMBCQ5 spectrum

2.1 SELHSQC and SELHSQCND spectra

Parameter sets: **awselhsqc** or **awselhsqcnd (+ getprosol)**

Pulse programmes: **awselhsqcgpsisp** or **awselhsqcndgpsisp**

Prior to running a SELHSQC experiment run a standard ^{13}C or DEPT experiment and determine the **O1** frequency in Hz of the ^{13}C signal to be selectively excited. Enter this value as **O2 (Hz)**.

TD = SI = 32 K.

SW = 12 ppm, O1P = 5 ppm. Adjust **SW** and **O1P** as required.

O2: frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, DS = 4 or 8.

D1 = 1 sec or other value of your choice.


D24 is automatically calculated from **CNST2** ($^1J_{\text{C-H}}$).

CNST2 = $^1J_{\text{C-H}}$; typically 125 to 160 Hz for $\text{sp}^3\text{-sp}^2$ carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have $^1J = 200\text{-}220$ Hz.

Check that gradients and shaped pulses are OK, including a **40000 usec p36:sp26** band selective **Q3.1000** pulse.

Set **receiver gain** using **RGA** (*Important!*).

Process with: **EFP** (applies **LB = 0.1-0.3 Hz**)



AVII-600 SELHSQC (lower) and SELHSQCND (upper) spectra determined for quinine in D₆-DMSO with selective excitation of the ^{13}C signal at 102.4 ppm (**O2 = 15534 Hz**).

2.2 SELHSQC-DIPS12 and SELHSQCND-DIPS12 spectra

Parameter sets: **awselhsqc-dipsi2** or **awselhsqcnd-dipsi2 (+ getprosol)**

Pulse programmes: **awselhsqcgpdigpsisp** or **awselhsqcgpdigpndsisp**

Prior to running a **SELHSQC-DIPS12** experiment run a standard ^{13}C or **DEPT** experiment and determine the **O1** frequency of the ^{13}C signal in **Hz** to be selectively excited. Enter this value as **O2 (Hz)**.

TD = **SI** = 32 K.

SW = 12 ppm, **O1P** = 5 ppm. Adjust **SW** and **O1P** as required.

O2: frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, **DS** = 4 or 8.

or **NS** x **TD0** scans where **TD0** = any positive number.

D1 = 1 sec or other value of your choice.

D9 = 80 msec or other time of your choice (6-160 msec).


D24 is automatically calculated from **CNST2** ($^1J_{\text{C-H}}$).

CNST2 = $^1J_{\text{C-H}}$; typically 125 to 160 Hz for $\text{sp}^3\text{-sp}^2$ carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have $^1J = 200\text{-}220$ Hz.

Check that gradients and shaped pulses are OK, including a **40000 usec p36:sp26** band selective **Q3.1000** pulse.

Set **receiver gain** using **RGA** (*Important!*).

Process with: **EFP** (applies **LB** = 0.1-0.3 Hz)



Expansion of the 7.0-8.5 ppm region of the AVII-600 SELHSQC-DIPS12 (*lower*) and SELHSQCND-DIPS12 (*upper*) spectra determined for quinine in $\text{D}_6\text{-DMSO}$ with selective excitation of the ^{13}C signal at 102.9 ppm (**O2** = 15534 Hz).

Correlated ^1H NMR signals observed in coupled SELHSQCND-DIPS12 spectra show ^1J , ^2J , or ^nJ $^{13}\text{C}\text{-}^1\text{H}$ couplings depending on the number of bonds between the selectively excited ^{13}C signal and correlated proton signals.

2.3 SELHMBC spectrum

Parameter set: (+ getprosol)

Pulse programme:

Prior to running a SELHMBC experiment run a standard ^{13}C or DEPT experiment and determine the **O1** frequency in Hz of the ^{13}C signal to be selectively excited. Enter this value as **O2 (Hz)**.

TD = SI = 32 K.

SW = 10 or 12 ppm, O1P = 5 or 6 ppm. Adjust **SW** and **O1P** as required.

O2: frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, DS = 4 or 8.

D1 = 1 sec or other value of your choice.


CNST2 = $^1J_{\text{C-H}}$; typically 125 to 160 Hz for $\text{sp}^3\text{-sp}^2$ carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have $^1J = 200\text{-}220$ Hz.

CNST13 = 8 Hz or other value of your choice.

D6 is autocalculated from **CNST13**

Check that gradients and shaped pulses are OK, including a **40000 usec p36:sp26** band selective **Q3.1000** pulse.

Process with **EFP** (applies **LB = 0.3 Hz**) *and* **MC**.



AVII-600 SELHMBCQ5 spectrum determined for the signal which occurs at 156.4 ppm.

If no signal is observed (due to T_2 or other relaxation issues) try the SELHMBCQ5.2 experiment which uses a **20000 usec Q5.1000** pulse. The SELHMBCQ5.2 experiment has a somewhat wider excitation window than that of the SELHMBCQ5 experiment which uses a **40000 usec Q5.1000** pulse