



KJM 9250

XHDEPT, HETCOR and COLOC Experiments on the AVII-600

Version 7.3

Topspin 3.2 Windows 7 AVII600



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2020

XHDEPT, HETCOR and COLOC Experiments on the AVII-600

1.0 Introduction

¹³C detected aw coded **XHDEPT, HETCOR and COLOC** parameter sets are set up with 2K acquired ¹³C points and 128 to 256 ¹H increments. ¹H and ¹³C spectral windows and their midpoints should be determined before setting up these experiments.

1.1 Processing

HXDEPT45, HXDEPT90 and HXDEPT135 experiments are phase sensitive experiments. These spectra should be phased **before** using the **abs1** and **abs2** commands.

The **HXDEPTQF, HETCOR and COLOC** experiments are absolute value experiments. Phasing is not required.

2.0 Experiments and Parameter Sets

The following ¹H detected **HXDEPT, HETCOR and COLOC** experiments and parameter sets have been set up on the AVI- 600 spectrometer.

- 2.1 HXDEPT45, HXDEPT90 and HXDEPT135** spectra
- 2.2 HXDEPTQF** spectrum
- 2.3 HETCOR (= HXCOQF)** spectrum
- 2.4 COLOC** spectrum

2.1 HXDEPT45, HXDEPT90 and HXDEPT135 spectra

Parameter sets: **awhxdept45**, **awhxdept90** or **awhxdept135 (+ getprosol)**
pulse programmes: **awhxdept45**, **awhxdept90** or **awhxdept135**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1P** = ¹³C spectral window midpoint in ppm.

Enter **O2P** = ¹H spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 128-256 (your choice).

NS = 8, 16 (multiple of 4 or 8 recommended), **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.


CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg 125-160 Hz).

p0 = **p3* 0.5, 1.0 or 1.5** for **45°, 90° or 135°** pulses.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K or 2K**
WDW(F1) = WDW(F2) = QSINE
SSB(F2) = SSB(F1) = 2
xfb, abs1 and abs2



AVII-600 HXDEPT135 spectrum of quinine in D₆-DMSO.

2.2 HXDEPTQF

Parameter set: **awhxdeptqf (+ getprosol)**
pulse programmes: **hxdeptqf**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1P** = ¹³C spectral window midpoint in ppm.

Enter **O2P** = ¹H spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 128-256 (your choice).

NS = 8, 16 (multiple of 4 or 8 recommended), **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.


CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg 125-160 Hz).

p0 = 45 degree pulse for a DEPT45-like outcome.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K or 2K**
WDW(F1) = WDW(F2) = SINE
SSB(F2) = SSB(F1) = 0
xfb, abs1 and abs2



AVII-600 HXDEPTQF spectrum of quinine in D₆-DMSO.

2.3 HETCOR spectrum

Parameter set: **awhetcor (+ getprosol)**

Pulse programme: **hxcoqf**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1P** = ¹³C spectral window midpoint in ppm.

Enter **O2P** = ¹H spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 128-256 (your choice).

NS = 8, 16 (multiple of 4 or 8 recommended), **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg: 125-160Hz).

CNST11 = **3** (used to auto calculate **D3**).

Type **ased** (enter) and review parameters used in the job.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K or 2K**

WDW(F1) = WDW(F2) = SINE

SSB(F2) = SSB(F1) = 0

xfb, abs1 and abs2



AVII-600 **HETCOR** spectrum of quinine in D₆-DMSO.

2.4 COLOC spectrum

Parameter sets: **awcoloc (+ getprosol)**

Pulse programme: **colocqf**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1P** = ¹³C spectral window midpoint in ppm.

Enter **O2P** = ¹H spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 128-256 (your choice).

NS = 8, 16 (multiple of 4 or 8 recommended), **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

D6 = 50-80 msec or other delay for ⁿJ evolution.

D18 = 30 msec or other delay for ⁿJ evolution.

Parameter set **D6** and **D18** delay values are those for ⁿJ = ~10 Hz.

Type **ased** (enter) and review parameters used in the job.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K or 2K**

WDW(F1) = WDW(F2) = SINE

SSB(F2) = SSB(F1) = 0

xfb, abs1 and abs2



AVII-600 COLOC spectrum of quinine in D₆-DMSO optimized for ⁿJ = ~10 Hz

