

Safe Operation Procedures – FTIR

Name of department: Department of Chemistry (Kjemisk Institutt)

Date: 10/3-2015, Per-Anders Hansen

Only trained personnel are allowed to use the FTIR! Contact Per-Anders Hansen for training: p.a.hansen@kjemi.uio.no.

Sample positioning: Make sure you don't put any heavy weights on the platform when positioning the sample. Anything over 500 g is considered heavy. This includes leaning on it!

- **If in doubt about anything, contact the person responsible for the photospectrometer!**